fr	NOTICE D'UTILISATION – PRECIVIA®	2
en	INSTRUCTION FOR USE – PRECIVIA®	3
de	GEBRAUCHSANWEISUNG – PRECIVIA®	4
it	ISTRUZIONI PER L'USO – PRECIVIA®	6
es	INSTRUCCIONES DE USO – PRECIVIA®	7
pt	INSTRUÇÕES DE USO – PRECIVIA®	8
nl	GEbruIKSAANWIJZING – PRECIVIA®	10
da	BRUGSVEJLEDNING – PRECIVIA®	11
fi	KÄYTTÖOHJEET – PRECIVIA®	12
sv	BRUKSANVISNING – PRECIVIA®	13
el	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ – PRECIVIA®	15
et	KASUTUSJUHEND – PRECIVIA®	16
lv	LIETOSANAS INSTRUKCIJA – PRECIVIA®	17
lt	NAUDOJIMO INSTRUKCIJOS – PRECIVIA®	19
pl	INSTRUKCJA UŻYTKOWANIA – PRECIVIA®	20
sl	NAVODILA ZA UPORABO – PRECIVIA®	21
sk	NAVOD NA POUZIVANIE – PRECIVIA®	23
cs	NAVOD K POUZITI – PRECIVIA®	24
ro	INSTRUCȚIUNI DE UTILIZARE – PRECIVIA®	25
bg	ИНСТРУКЦИИ ЗА УПОТРЕБА – PRECIVIA®	27
hr	UPUTE ZA UPOTREBU – PRECIVIA®	28
hu	HASZNALATI UTASITAS – PRECIVIA®	30
no	BRUKSANVISNING – PRECIVIA®	31
tr	KULLANIM – PRECIVIA®	32
Zh-tw	使用說明 – PRECIVIA®	33
Illustrations		35
Symbols.....		36
IFU-2138 367 11 C		

FCI S.A.S.
20-22, rue Louis Armand
75015 PARIS – FRANCE

 : + 33 (0) 1 53 98 98 98

 : + 33 (0) 1 53 98 98 99

@ : fci@fci.fr

 : www.fci.fr

Distributed in U.S.A. by:

FCI Ophthalmics Inc.
30 Corporate Park Drive
Suite # 310/320
Pembroke, MA 02359

 : 800-932-4202

 : 781-826-9062

@ : info@fci-ophthalmics.com

 : www.fci-ophthalmics.com

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Product scope of this IFU:

- S5.200.I: PRECIVIA® I
- S5.200.P: PRECIVIA® P

PRECIVIA® is a medical device made of transparent polycarbonate. The device is composed of:

- the body of the device,
- the flange in contact with the eyeball
- the tube guide to ensure a fixed injection angle of 28°, a fixed distance at 3.5 mm from the limbus, and a secure injection depth (between 5 and 15 mm, the actual depth depending on the needle used).

The device is sterilized using ethylene oxide in its final simple pouch packaging. A green indicator shows that the product followed a sterilization cycle validated by FCI. The device is supplied by box of 25. The device is for single use only and is not resterilizable.

Intended purpose

PRECIVIA® is a device that standardizes, secures and simplifies intravitreal injections and makes them reproducible, by guiding the orientation and depth of the needle.

Indications

Medical conditions: Standardizes intravitreal injections for (non-exhaustive list):

- AMD (neovascular age related macular degeneration)
- CSME (clinically significant macular edema)
- PDR (proliferative diabetic retinopathy)
- Retinal Vein Occlusions
- Uveitis
- CME (cystoid macular edema)

Intended patient population: adults

Body parts in contact with the device: eyes, cornea, under upper and lower eyelids

Intended users: ophthalmologists or trained nurses only

Intended use environment: operating room

Contraindications

As PRECIVIA® is set under the lid, its use may be avoided in patients who have undergone previous filtering glaucoma surgery to prevent bleeding-related endophthalmitis.

Warnings

The choice of the PRECIVIA® model to use must be made from the indications in the paragraph "Compatibility with other devices": a check of the compatibility should be performed by the user before its use on the patient. Using an incompatible needle poses a risk to the patient regarding the depth of injection.

The placement of PRECIVIA® may be difficult and painful if patients tighten their lids excessively during the insertion.

The guaranteed maximum residual rate of ethylene oxide the device can release to the patient is 60 µg per device. The residual quantity of ethylene oxide corresponds to the maximum quantity of residue that may leach over the lifetime of use.

The device is for single use and must not be resterilized. Reuse and/or resterilization may compromise device performance, which could cause serious harm to the patient's health and safety.

Check packaging integrity before use, and do not use the device in the event of packaging being damaged, unintentionally opened before use, or if it was exposed to storage or transport conditions outside of those specified.

Do not use the product if the indicator is not green, as it could endanger the patient. Do not use after the expiration date shown on the package.

In case the user identifies faulty labelling information, it shall be relayed to FCI S.A.S. and the device should not be used if the patient's safety is threatened.

Precautions for use

The device must be stored at temperatures between 5°C/45°F and 35°C/95°F and humidity between 30% and 75%, away from sunlight.

The device may be exposed to temperature variations between -15°C/5°F and 50°C/122°F and humidity variations between 15% and 90% during transport for 10 days.

Remove the device from its packaging in aseptic conditions.

Local anesthetic must be performed prior to the setting of PRECIVIA® and the injection of drugs. Topical antibiotics and ophthalmic viscoelastic device can be instilled to the eye following injection.

Adverse effects

As in any type of surgery, there are risks linked to the material or the evolution of the initial pathology. Potential complications associated with the use of the device include,

but are not limited to the following:

- subconjunctival hemorrhages
- corneal abrasion
- infection

Other risks that can occur are directly related to intravitreal injection and not the use of PRECIVIA®. These risks include (non-exhaustive list):

- Pain
- Retinal tear / detachment
- Cataract
- Loss of vision
- Increased IOP

Any complication that occurs in relation with the device should be notified to FCI S.A.S. and the competent authority of the Member State in which the user and/or the patient is established, or, for users and/or patients not subject to EU requirements, to FCI S.A.S. and to the applicable regulatory authority in conformity with the requirements of that regulatory authority.

Instructions for use

1. Please select the PRECIVIA® model suitable for the needle selected for the injection: the compatibilities between PRECIVIA® and the injection equipment are described below in the section "Compatibility with other devices".
2. Place PRECIVIA® on the patient's eye with the central orientation mark correctly lined with the limbus.
3. Proceed with a smooth downward pressure with a slight rotation ensuring a certain amount of conjunctival displacement and eye fixation.
4. Ask the patient to look in a defined area (the patient is then unable to see the needle).
5. Insert the needle into the eye using the guide tube of the PRECIVIA® device and inject the drug into the eye.
6. Remove the needle and the PRECIVIA® device using a reversed rotation

Performances, mode of action and clinical benefits

There is no clinical performance associated with the use of PRECIVIA®. The device's mode of action is mechanical guidance of the needle at a fixed angle (28°), a position (3.5 mm from the limbus) and a secure injection depth (between 5 mm to 15 mm). The benefit for the patient is that the intravitreal injection is simplified, safe and standardized.

Compatibility with other devices

PRECIVIA® P (S5.200.P) is intended to be used with syringes filled with solutions and needles 30 Ga x ½ inches, connected by a luer-lock connection.

PRECIVIA® I (S5.200.I) is compatible with insulin syringes with crimped needle

Information to be communicated to the patient

Any serious incident that occurs in relation with the device should be notified to FCI S.A.S. and the competent authority of the State in which the patient is established.

The user must inform the patient of the contraindications and potential complications related to the device that are set out in the present document

Waste disposal

With household waste or recycle	With hospital waste (contaminated products)
Outer box, instructions for use	PRECIVIA®
Sachet	

The disposal of this device does not present any physical risk.

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